

BiyoClamp[®]

MALE INCONTINENCE CLAMP



The most economical, comfortable and practical way to prevent male pattern urinary incontinence...



www.biyoclamp.com

Why BiyoClamp?



MALE INCONTINENCE

BiyoClamp is a clamp that controls involuntary urine leakage for men with urinary incontinence.

Thanks to the triangle at the bottom of the BiyoClamp, it controls the flow of urine by pressing on the urinary tract.




PRP, P-SHOT

BiyoClamp supports the treatment by being attached to the penis during the P-Shot procedure.



PELVIC FLOOR EXERCISES

BiyoClamp contributes to the treatment process by supporting pelvic floor exercises.



RADIOTHERAPY

It is attached to the penis during radiotherapy treatment to keep bladder full.

ADVANTAGES



PERSONAL HYGIENE

BiyoClamp does not cause odor and diaper rash like diapers.



RISK OF INFECTION

BiyoClamp leaves the urinary tract open and the risk of infection is almost non-existent, indwelling urinary catheters are attached directly to the urinary tract, so the risk of infection is quite high.

BiyoClamp does not affect activities such as sexual intercourse and exercise and provides the opportunity to insert and remove it when necessary. Permanent urinary catheters do not allow for these activities.



COMFORTABLE LIFE

BiyoClamp, provides a comfortable life all day thanks to its durable and lightweight design.



SUSTAINABLE AND ENVIRONMENTALLY FRIENDLY

It is environmentally friendly and does not infect household.

SELECTING THE CORRECT SIZE



How to Choose the Suitable Size?

BioClamp product measurements are designed according to the length (thickness) of the penis circumference. In order to choose the right product size, it is necessary to measure the thickness of the area between the glans penis and the groin with a tape measure

Choosing the appropriate size is important!

Choosing the appropriate size is important before using BiyoClamp! The use of inappropriately large sizes may be insufficient to control urine flow. Using an inappropriately small size may cause both irritation in the penis and disruption of blood flow.



SIZES

SMALL :8-9-10cm
MEDIUM :11-12cm
LARGE :13cm +



HOW TO USE?

Choosing the appropriate size is important before using BiyoClamp!
The use of inappropriately large sizes may be insufficient to control urine flow. Using an inappropriately small size may cause both irritation in the penis and disruption of blood flow.
BiyoClamp should always be applied to a clean and dry penis. In this way, it can be used for at least 6 months.



The urinary tract under the penis should be placed on the penis in a way that corresponds to the triangle at the bottom of the BiyoClamp.
BiyoClamp has a 5-stage locking system and should always be used at the level where urine flow stops.
Urine should be emptied by opening the BiyoClamp every 2 hours at the latest and then a break should be taken for 15-20 minutes without wearing it.



BiyoClamp[®]

MALE INCONTINENCE CLAMP

For further information,

Please contact: info@biyotip.net

The specifications shown in this leaflet are for information only and are not, under any circumstances, of a contractual nature.

Communication

Yıldırım Beyazıt District Aşık Veysel Boulevard

Erciyes Technopark 3.Block No:67/20 38010

Melikgazi Kayseri

+90.352 225 24 95

+90 533 720 67 73

info@biyotip.net

BIYOTIP

TIBBİ ÜRÜNLER SAN.VE TİC.LTD.ŞTİ